As an accredited laboratory, this laboratory is entitled to use the following accreditation symbol.





## Schedule of Accreditation

## Accreditation Scheme for Medical/Clinical Laboratories Sri Lanka Accreditation Board for Conformity Assessment

## Accreditation Number: ML 019-01

## Ohlums Clinic & Laboratories (Pvt)Ltd No.132, Dr N M Perera Mawatha Cotta road, Borella Colombo-08

Scope of Accreditation: Performing Medical Tests under the fields of Clinical Biochemistry, Clinical Pathology, Haematology and Serology

The laboratory is accredited for the following tests

SI	Field of Testing	Test	Test Method	Test Instrument	Range of Testing	MU
No						(±)
01	Clinical	Creatinine	Kinetic	Chem Well 2910	0.176-17 mg/dl	Level 1: 8.82
	Biochemistry		modified Jaffe			Level 2: 9.53
		Blood Urea	U.V		2 -300 mg/dl	Level 1: 9.14
						Level 2: 6.97
		AST/SGOT	Kinetic U.V		2 -700 U/I	Level 1: 6.23
						Level 2: 5.96
		ALT/SGPT	Kinetic U.V		4 -600 U/I	Level 1: 9.14
						Level 2: 7.38
02	Clinical	HCG	Rapid	Strip	1.Detected	
	Pathology	(Urine)	Chromatograph	-	2.Non Detected	-

SI No	Field of Testing	Test	Test Method	Test Instrument	Analytical Range	MU (±)			
03	Haematology	Full Blood Count							
		Haemoglobin	Non-Cyanide Haemoglobin Analysis	Sysmex heamatalogy	0.1-25.0 g/dl	0.4			
		PCV (Packed cell volume)	Calculated from RBC & HGB	Analyzer Kx-21	10.0-60.0 (%)	2.6			
		RBC(Red blood cell)	Direct current resistance changes		$0.3 - 7.0 (10^{6}/\mu L)$	0.18			
		MCV(Mean corpuscular volume)	Calculated from RBC & HCT		N/A	3.9			
		MCH(Mean corpuscular Haemoglobin)	Calculated from RBC & HGB		N/A	1.6			
		MCHC(Mean corpuscular Haemoglobin concentration)	Calculated from HCT & HGB		N/A	2.6			
		Platelet count	Direct current resistance changes	-	10-999 (10 <sup>6</sup> /μL)	33			
		WBC count(White blood cell count)			1.0 –999 (10 <sup>6</sup> /μL)	0.5			
		Differential	Manual	Microscope	0-100	-			
		ESR	Westegren	Westegren Tube & E.S.R Rack	0-180	-			
04	Serology	HIV 1& 2	ELISA	EVOLIS TWIN	Positive				
		Hepatitis B Hepatitis C		PLUS Automated	Negative Indeterminate	-			
		VDRL	Carbon Antigen	Rotary Shaker White Slide (with separated square) & Micro pipette	Reactive Non-reactive Weakly reactive	-			

Director /CEO Sri Lanka Accreditation Board for Conformity Assessment