



Valid from 04 June 2025
to 03 June 2029
Issued on 23 June 2025

As an accredited laboratory, this laboratory is entitled to
use the following accreditation symbol.



ISO 15189
ML 051-01

Schedule of Accreditation

Accreditation Scheme for Medical/Clinical Laboratories
Sri Lanka Accreditation Board for Conformity Assessment

Accreditation Number: ML 051-01

Histopathology Laboratory
Lady Ridgeway Hospital for Children
Dr. Danister De Silva Mawatha
Colombo 08.

Scope of Accreditation: Performing Medical/Clinical testing under the fields of
Histopathology

The laboratory is accredited for the following tests.

Sl	Field of Testing	Test	Test Method	Test Instrument
1	Histopathology	Tissue processing and staining with H & E (Tissue Sample)	LRH/H/M/04/2024/Tissue process manual and LRH/H/SOP/11/SOP for H & E	Cut Up-Grossing station Processing Vacuum Processor Embedding-Embedding console Cutting Microtome - Done Manually Slide Warmer-Done Manually Staining - Automated Slide Stainer Automated Cover Slipper-Automated
2		Tissue processing and Special stain PAS (periodic acid Schiff) (Tissue sample)	LRH/H/M/04/2024/Tissue process manual and LRH/H/SOP/13/SOP for PASD	Cut Up-Grossing station Processing Vacuum Processor Embedding-Embedding console Cutting Microtome - Done Manually Slide Warmer-Done Manually Staining with PAS

Sl	Field of Testing	Test	Test Method	Test Instrument
3	Histopathology	Tissue processing and Special stain PASD (periodic acid shiff diastase) (Tissue Sample)	LRH/H/M/04/2024/Tissue process manual and LRH/H/SOP/13/SOP for PASD	Cut Up-Grossing station Processing Vacuum Processor Embedding-Embedding console Cutting Microtome - Done Manually Slide Warmer-Done Manually Staining with PAS/PASD
4		Tissue processing and Special stain-Masson's trichrome (Tissue sample)	LRH/H/M/04/2024/Tissue process manual and LRH/H/SOP/14/SOP for Massons	Cut Up-Grossing station Processing Vacuum Processor Embedding-Embedding console Cutting Microtome - Done Manually Slide Warmer-Done Manually Staining with Massons
5		Tissue processing and Special stain- silver stain (Tissue sample)	LRH/H/M/04/2024/Tissue process manual and LRH/H/SOP/30/SOP for silver staining	Cut Up-Grossing station Processing Vacuum Processor Embedding-Embedding console Cutting Microtome - Done Manually Slide Warmer-Done Manually Staining with Silver
6		Tissue processing and Special stain- Alcian blue (Tissue sample)	LRH/H/M/04/2024/Tissue process manual and LRH/H/SOP/16/SOP for Alcian blue	Cut Up-Grossing station Processing Vacuum Processor Embedding-Embedding console Cutting Microtome - Done Manually Slide Warmer-Done Manually Staining with Alcian blue
7		Tissue processing and Special stain- Ziehl Neelsen (Tissue sample)	LRH/H/M/04/2024/Tissue process manual and LRH/H/SOP/15/SOP for Ziehl Neelsen	Cut Up-Grossing station Processing Vacuum Processor Embedding-Embedding console Cutting Microtome - Done Manually Slide Warmer-Done Manually appropriate Staining
8		Tissue processing and Special stain- Wade-Fite (Tissue sample)	LRH/H/M/04/2024/Tissue process manual and LRH/H/SOP/17/SOP for Wade-Fite	Cut Up-Grossing station Processing Vacuum Processor Embedding-Embedding console Cutting Microtome - Done Manually Slide Warmer-Done Manually appropriate Staining

Sl	Field of Testing	Test	Test Method	Test Instrument
9	Histopathology	Cytology stain with H & E (Fluid sample)	LRH/H/SOP/04/SOP for routine cytology sample collection, transport and processing LRH/H/M/05/2024/FNA Manual	Cytocentrifuge/Centrifuge Fixing Stain
10		FNA (Fine needle aspiration test) (Fluid sample)	LRH/H/M/05/2024/FNA Manual	Manually aspirate Fixed and Stain

Acting Director/CEO
Sri Lanka Accreditation Board for Conformity Assessment